



#### UV Fluorescence & Colorimetry

# Essential V6

The Essential V6 imaging platform is Uvitec's most compact, yet high-end gel documentation system for DNA, RNA and protein gel imaging, boasting a 3-megapixel, 16-bit sensor offering massive 65,535 gray scale ideal for lab imaging routine, research-level data collection and accurate protein quantitative analysis

**Dims** > Height: 570 mm - Width: 390 mm - Depth: 475 mm **Weight** > With Transilluminator: 36 Kg - Without: 30 Kg

## Practicality

- ▶ 1-click to image > effortless acquisition, in no time
- ▶ No manual control > hands-off, automated routines
- ▶ Fully slide-out tray > hassle-free sample positioning
- ▶ **All-in-one** > acquisition editing analysis, straight away

## **Applications**

- ▶ **Uvipure© technology** > enhanced UV for EtBr and all safe stains
- ▶ **Multi-functions** > optional screens for extended applications
- Quantitative analysis > against one or several band(s) of reference
- ► **Fully-customizable** > 100% adjustable parameters and protocols

### Optics

- ▶ **3,000,000 pixels** > superb, publication-level images
- ▶ 65,535 gray levels > research-level protein quantification
- **6x zoom lens** > responsive optics for a versatile imaging area
- ▶ Passive cooling > reduced noise for clearer images

#### Darkroom

- ▶ Stainless steel > long-lasting, robust instrument
- ▶ **Epoxy paint** > chemical-resistant darkroom
- ▶ Compact design > small footprints and easy fit
- ▶ **Lightproof** > no parasite signals, no light leaks

► Fluorescence

DNA and RNA gels with fluorescent stains

Visible Imaging

Optional Colorimetry and Protein gels